



VERTICAL AUTOCLAVE



YXQ-LB-SII Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

VERTICAL AUTOCLAVE

YXQ-LB-SII Series

Product Features

- Made of SUS304
- Hand wheel type quick-open cover;
- Microcomputer intelligent control;
- Self-expanding type seal ring;
- Exhaust cold air automatically;
- Self-control when lack of water.
- Self-exhaust when over pressure;
- Pressure safety interlock;
- Over temperature protection;
- Buzzer reminds when sterilization ends;
- Equipped with G1/2A internal thread sample test hole.



Model	YXQ-LB-30SII	YXQ-LB-50SII	YXQ-LB-75SII	YXQ-LB-100SII
Volume [L]	29	50	75	100
Power [kW]	2.3	3.2	3.5	3.5
Power Supply [V, Hz]	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%
Max. Working / Design Temp. [°C]	135 / 139	135 / 139	135 / 139	135 / 139
Max. Working / Design Pressure [MPa]	0.22 / 0.25	0.22 / 0.25	0.22 / 0.25	0.22 / 0.25
Timer Range [min]	4~120	4~120	4~120	4~120
Internal Dimensions [mm]	Ø300×485	Ø400×630	Ø400×720	Ø400×830
Housing Dimensions W×D×H [mm]	445×445×900	485×485×1120	485×485×1190	485×485×1290
Basket Dimensions [mm]	Ø280×340	Ø360×230	Ø360×270	Ø360×230
Number of Baskets [pc]	1	2	2	3
Net Weight [kg]	54	68	71	75
Gross Weight [kg]	64	82	85	89
Package Size [mm]	570×570×1060	610×610×1290	610×610×1360	610×610×1490